


Product Datasheet

p70 S6 Kinase (phospho-Ser371) antibody (orb6620)

Description	p70 S6 Kinase (phospho-Ser371) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IHC-P 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO whole cell...