

Product Datasheet

Mast Cell Chymase Immunizing peptide (orb66255)

Description	Mast Cell Chymase Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb66255 (Mast Cell Chymase antibody)
Protein Sequence	QFRHPKYNTSTLHHD
Expiration Date	12 months from date of receipt.