

Product Datasheet

Methyl-PABP(Arg455/460) antibody (orb6630)

Description	Methyl-PABP(Arg455/460) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised metylpeptide derived from human PABP around the methylation site of Arg455/Arg460 (GAI(Me-R)PAAP(Me-R)PP)
Target	PABP (Methyl Arg455/460)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75 kDa
Uniprot ID	P11940
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.