

Product Datasheet

PAK4 (phospho-Ser474) antibody (orb6634)

Description	PAK4 (phospho-Ser474) antibody
Species/Host	Rabbit
Reactivity	Feline, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PAK6 around the phosphorylation site of Ser560 (RK(p-S)LV)
Target	Phospho-PAK6 (Ser560)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75 kDa
Uniprot ID	Q9NQU5
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.