



Product Datasheet

Histone H1oo (AcK163) antibody (orb665467)

Description	Rabbit polyclonal antibody to Histone H1oo (Acetyl-K163)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K163 of human Histone H1oo protein. The exact sequence is proprietary.
Target	H1FOO
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8IZA3 , Q8VIK3
Entrez	171506 , 132243
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H1oo (A...



Immunohistochemical analysis of Histone ...