

Product Datasheet

Histone H2A (AcK95) antibody (orb665469)

Description

Rabbit polyclonal antibody to Histone H2A (Acetyl-

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K95 of human Histone H2A protein. The exact sequence is proprietary.

Target

HIST1H2AG; HIST1H2AI; HIST1H2AK; HIST1H2AL; HIST1H2AM

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P02262](#), [P22752](#), [P0C0S8](#)

Entrez


[319170](#), [319172](#), [8329](#), [8330](#), [8332](#), [319164](#), [319166](#), [8336](#), [319167](#), [319191](#), [8969](#), [319165](#), [665433](#), [319171](#)

Dilution Range


WB: 1-500-1000

Expiration Date

12 months from date of receipt.



Western blot analysis of Histone H2A (Acetylated)



Direct ELISA antibody dose-response curve