

Product Datasheet

Histone H2B (AcK85) antibody (orb665470)

Description

Rabbit polyclonal antibody to Histone H2B (Acetyl-



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot analysis of Histone H2B (Acetylated)

Immunogen

KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K85 of human Histone H2B protein. The exact sequence is proprietary.

Target

HIST1H2BC; HIST1H2BE; HIST1H2BF; HIST1H2BG; HIST1H2BI; HIST1H2BD; HIST1H2BH; HIST1H2BK; HIST1H2BL; HIST1H2BM; HIST1H2BN; HIST2H2BF; H2BFS

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P62807](#), [P23527](#), [Q99879](#), [Q6ZWH9](#), [P58876](#), [P06899](#), [Q99877](#), [Q64475](#), [P33778](#), [Q93079](#), [Q99880](#), [Q64478](#), [O60814](#)

Entrez

[8344](#), [8346](#), [85236](#), [319178](#), [319182](#), [319179](#), [8340](#), [8339](#), [8347](#), [8342](#), [440689](#), [3017](#), [8343](#), [8345](#), [8341](#)

Dilution Range

WB: 1-500-1000

Expiration Date

12 months from date of receipt.