



Product Datasheet

INSL4 antibody (orb665585)

Description	Rabbit polyclonal antibody to INSL4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human INSL4. The exact sequence is proprietary.
Target	INSL4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q14641
Entrez	3641
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200, IF: 1-100-500, IP: 1-10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of INSL4 expression...



Immunohistochemical analysis of INSL4 staining...



Immunofluorescent analysis of INSL4 staining...