



Product Datasheet

SPTBN5 antibody (orb665609)


Description	Rabbit polyclonal antibody to SPTBN5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human SPTBN5. The exact sequence is proprietary.
Target	SPTBN5
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9NRC6
Entrez	51332
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200, ICC: 1-100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of SPTBN5 expressi...



Immunohistochemical analysis of SPTBN5 S...



Immunocytochemistry analysis of SPTBN5 S...