


Product Datasheet

CHRM2 antibody (orb665625)

Description	Rabbit polyclonal antibody to CHRM2
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CHRM2. The exact sequence is proprietary.
Target	CHRM2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P08172, P10980, Q9ERZ4
Entrez	1129, 243764, 81645, 102551851
Dilution Range	WB: 1-500-1000, IF: 1-100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of CHRM2 expression...



Immunofluorescent analysis of CHRM2 staining...