

Product Datasheet

Synaptotagmin (pS309/306) antibody (orb665649)

Description

Rabbit polyclonal antibody to Synaptotagmin

Species/Host

Rabbit

Reactivity

Human, Mouse, Primate, Rat

Conjugation

Unconjugated

Tested

IH, WB

Applications

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S309/306 of human Synaptotagmin protein. The exact sequence is proprietary.

Target

SYT1; SYT2

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P46097](#), [P21579](#), [P29101](#), [P21707](#), [P46096](#), [Q8N9I0](#)

Entrez


[20980](#), [20979](#), [6857](#), [24805](#), [127833](#), [25716](#)

Dilution Range


WB: 1-500-1000, IHC-P: 1-50-200

Expiration Date

12 months from date of receipt.



Western blot analysis of Synaptotagmin (...)



Immunohistochemical analysis of Synaptot...