

Product Datasheet

Alpha-tubulin (pY272) antibody (orb665690)

Description

Rabbit polyclonal antibody to Alpha-tubulin

Species/Host

Rabbit

Reactivity

Human, Mouse, Primate, Rat, Zebrafish

Conjugation

Unconjugated

Tested

IH, WB

Applications

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y272 of human Alpha-tubulin protein. The exact sequence is proprietary.

Target

TUBA1A

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P68369, Q71U36, P68366, P68370, Q13748, Q9BQE3, Q6AYZ1, Q6P9V9, P05213, P68373, P68363, Q6PEY2

Entrez


113457, 22146, 22142, 500929, 64158, 7278, 100909441, 22143, 7277, 112714, 84790, 10376, 7846

Dilution Range


WB: 1-500-1000

Expiration Date

12 months from date of receipt.



Western blot analysis of Alpha-tubulin (...)



Immunohistochemical analysis of Alpha-tu...