



Product Datasheet

SOX8 antibody (orb665742)

Description	Rabbit polyclonal antibody to SOX8
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human SOX8. The exact sequence is proprietary.
Target	SOX8; SOX9; SOX17; SOX18
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P48436 , Q9H6I2 , Q61473 , Q04886 , Q04887 , P35713 , P43680 , P57073
Entrez	20682 , 64321 , 20681 , 54345 , 30812 , 20672 , 6662 , 20671
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of SOX8 expression...



Immunohistochemical analysis of SOX8 sta...