



Product Datasheet

RhoBTB1 antibody (orb665757)

Description	Rabbit polyclonal antibody to RhoBTB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human RhoBTB1. The exact sequence is proprietary.
Target	RHOBTB1; RHOBTB2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9DAK3 , O94844 , Q9BYZ6 , Q91V93
Entrez	9886 , 246710 , 69288 , 23221
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200, IF: 1-100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of RhoBTB1 express...



Immunohistochemical analysis of RhoBTB1

...



Immunofluorescent analysis of RhoBTB1 st...