


Product Datasheet

PRKD1/2/3 antibody (orb665758)

Description	Rabbit polyclonal antibody to PRKD1/2/3
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PRKD1/2/3. The exact sequence is proprietary.
Target	PRKD1; PRKD2; PRKD3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9WTQ1 , Q8K1Y2 , Q15139 , Q5XIS9 , Q62101 , Q8BZ03 , Q9BZL6 , O94806
Entrez	85421 , 25865 , 18760 , 292658 , 75292 , 101540 , 5587 , 23683
Dilution Range	WB: 1-500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PRKD1/2/3 expression.