



## Product Datasheet

### Frizzled 5 antibody (orb665816)

|                            |                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Frizzled 5                                                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                        |
| <b>Tested Applications</b> | IF, WB                                                                                                                                              |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Frizzled 5. The exact sequence is proprietary.           |
| <b>Target</b>              | FZD5; FZD8                                                                                                                                          |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                           |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                           |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                    |
| <b>Note</b>                | For research use only                                                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                                                                          |
| <b>Source</b>              | Rabbit                                                                                                                                              |
| <b>Uniprot ID</b>          | <a href="#">Q498S8</a> , <a href="#">Q13467</a> , <a href="#">Q61091</a> , <a href="#">Q8CHL0</a> , <a href="#">Q9EQD0</a> , <a href="#">Q9H461</a> |
| <b>Entrez</b>              | <a href="#">100909849</a> , <a href="#">7855</a> , <a href="#">317674</a> , <a href="#">14370</a> , <a href="#">8325</a> , <a href="#">14367</a>    |
| <b>Dilution Range</b>      | WB: 1-500-1000, IF: 1-100-500                                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                     |



Western blot analysis of Frizzled 5 expression.



Immunofluorescent analysis of Frizzled 5 expression.