

Product Datasheet

Actin-pan (AcK52/50/51) antibody (orb665845)

Description

Rabbit polyclonal antibody to Actin-pan (Acetyl-

kDa A
130

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western
blot
analysis of
Actin-pan
(Acet...

Immunogen

KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K52/50/51 of human Actin-pan protein. The exact sequence is proprietary.

Target

ACTA2; ACTB; ACTC1; ACTG1; ACTG2; ACTA1

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P63267, P60710, P63260, P68032, P68033, P63259, P60709, P63261, P62736, P68035, P68133, P62737, P62738, P60711

Entrez

71, 11475, 11461, 70, 81633, 58, 11465, 29275, 60, 100361457, 11464, 81822, 59, 72

Dilution Range

WB: 1-500-1000

Expiration Date

12 months from date of receipt.

