



Product Datasheet

NMDAR1 (pS890) antibody (orb665849)

Description	Rabbit polyclonal antibody to NMDAR1 (Phospho-S890)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S890 of human NMDAR1 protein. The exact sequence is proprietary.
Target	GRIN1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q05586 , P35438
Entrez	2902 , 14810
Dilution Range	WB: 1-500-1000, IHC-P: 1-50-200, IF: 1-100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of NMDAR1 (Phospho-S890)



Immunohistochemical analysis of NMDAR1 (Phospho-S890)



Immunofluorescent analysis of NMDAR1 (Phospho-S890)