



Product Datasheet

MTMR9 antibody (orb665907)


Description	Rabbit polyclonal antibody to MTMR9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant protein of human MTMR9. The exact sequence is proprietary.
Target	MTMR9
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9Z2D0 , Q96QG7
Entrez	66036 , 210376
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MTMR9 expressio...



Immunohistochemical analysis of MTMR9 stai...



Immunofluorescent analysis of MTMR9 stai...