



Product Datasheet

IL-25 antibody (orb665917)


Description	Rabbit polyclonal antibody to IL-25
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant protein of human IL-25. The exact sequence is proprietary.
Target	IL25
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9H293
Entrez	64806
Dilution Range	WB: 1-500-1000, IHC-P: 1-100-200, IF: 1-100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of IL-25 expressio...



Immunohistochemical analysis of IL-25 st...



Immunofluorescent analysis of IL-25 stai...