



Product Datasheet

ACAA2 antibody (orb665953)


Description	Rabbit polyclonal antibody to ACAA2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human ACAA2. The exact sequence is proprietary.
Target	ACAA2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8BWT1 , P13437 , P42765
Entrez	10449 , 170465 , 52538
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of ACAA2 expression...



Immunohistochemical analysis of ACAA2 staining...



Immunofluorescent analysis of ACAA2 staining...