


Product Datasheet

E2 Tag antibody (orb66646)

Description	Mouse monoclonal to E2 Tag.E2 tag peptide GVSSTSSDFRDR conjugated to KLH. E2 Tag antibody can recogn...
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:3000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of E2 fusion prote...