

Product Datasheet

PDGFR beta (phospho-Tyr1021) antibody (orb6665)

Description	PDGFR beta (phospho-Tyr1021) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDGF Receptor beta around the phosphorylation site of Tyr1021 (ND(p-Y)II)
Target	Phospho-PDGFR Receptor beta (Tyr1021)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190 kDa
Uniprot ID	P09619
Expiration Date	12 months from date of receipt.