



Product Datasheet

mCherry antibody (orb66657)

Description	mCherry antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Human mCherry recombinant protein (orb178365) was used as immunogen.
Target	mCherry
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	At -20°C for one year. After reconstitution at 4°C for two weeks. It can also be aliquotted and stored frozen at -20°C for a longer time.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Monoclonal antibodies are purified by peptide affinity chromatography
MW	34 kDa
Dilution Range	WB: 1 µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of protein 20 ng (lane 1), p...



Western blot analysis of mcherry protein...