


Product Datasheet

beta Tubulin antibody (orb66666)


Description	Mouse monoclonal to beta Tubulin. Microtubules are constituent parts of the mitotic apparatus, cilia,...
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Monkey, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:10000-500000, IHC-P: 1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 (1) , Rat ...



Mouse brain (3)



Chicken lung (4) and Rabbit testis (5)