


Product Datasheet

IL17A antibody (orb66667)

Description	Mouse monoclonal to IL17A. IL-17A is a proinflammatory cytokine produced by activated T cells. IL-17A...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of human IL-17A
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	17.5
Uniprot ID	Q16552
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...