


Product Datasheet

CD31 antibody (orb66672)

Description	Mouse monoclonal antibody to CD31
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse CD31
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	83
Uniprot ID	Q08481
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...