


Product Datasheet

GR1 antibody (orb66674)

Description	Mouse monoclonal to GR1. Granulocyte-differentiation antigen-1 (GR-1) is also known as Lymphocyte ant...
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse GR1
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	12
Uniprot ID	P35461
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...