

Product Datasheet

IL20 antibody (orb66681)

Description	Mouse monoclonal to IL20. Interleukin 20 (IL20) is a cytokine that signals through signal transducer ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	C-terminal region of human IL20
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	20
Uniprot ID	Q9NYY1
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.