


Product Datasheet

CNPASE antibody (orb66684)

Description	Mouse monoclonal to CNPASE.2', 3'-cyclic-nucleotide 3'-phosphodiesterase (CNPASE) belongs to the cyc...
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse CNPASE
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	47
Uniprot ID	P16330
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...