


Product Datasheet

GST antibody (orb66695)

Description	Mouse monoclonal to GST.Glutathione S-transferase class-mu 26 kDa isozyme (GST26) is an enzyme that ...
Species/Host	Mouse
Reactivity	S. japonicum
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	immunizing animals with the full-length GST26 from Schistosoma japonicum (Blood fluke)
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	26
Uniprot ID	P08515
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Schistosoma japonicum.