


Product Datasheet

PRMT8 antibody (orb66696)

Description	Rabbit polyclonal to PRMT8. Glutathione S-transferase class-mu 26 kDa isozyme (GST26) is an enzyme th...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PRMT8
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	Q9NR22
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...