


Product Datasheet

CCL1 antibody (orb66705)

Description	Rabbit polyclonal to CCL1.C-C motif chemokine 1 (CCL1) belongs to the intercrine beta (chemokine CC)...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse CCL1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	10
Uniprot ID	P10146
Dilution Range	WB:1:1000-1:5000, ELISA:Greater than 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...