


Product Datasheet

SATB1 antibody (orb66710)

Description	Rabbit polyclonal to SATB1. Special AT-rich sequence-binding protein 1 (SATB1) is a chromatin organiz...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SATB1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	86
Uniprot ID	Q01826
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...