


Product Datasheet

DNMT3A antibody (orb66719)

Description	Rabbit polyclonal to DNMT3A. DNA (cytosine-5)-methyltransferase 3A (DNMT3A) is required for genome wi...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human DNMT3A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	102
Uniprot ID	Q9Y6K1
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...