


Product Datasheet

DNMT3A antibody pair (orb66721)

Description	Rabbit polyclonal to DNMT3A. 50 μl DNMT3A (N-term2) &50 μl DNMT3A (N-term). vPair antibodies r...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human DNMT3A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	102
Uniprot ID	Q9Y6K1
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...