


Product Datasheet

HDAC1 antibody (orb66723)

Description	Rabbit polyclonal to HDAC1. Histone deacetylase 1 (HDAC1) is responsible for the deacetylation of lys...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human HDAC1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	55
Uniprot ID	Q13547
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...