


Product Datasheet

DNMT1 antibody (orb66726)

Description	Rabbit polyclonal to DNMT1. DNA (cytosine-5)-methyltransferase 1 (DNMT1) plays a role in methylating ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human DNMT1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	183
Uniprot ID	P26358
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...