


Product Datasheet

HAT1 antibody (orb66733)

Description	Rabbit polyclonal to HAT1. Histone acetyltransferase 1 (HAT1) is involved in the rapid acetylation of...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human HAT1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	O14929
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...