


Product Datasheet

MUC20 antibody (orb66735)

Description	Rabbit polyclonal to MUC20. Mucin 20 (MUC20) is a member of the mucin protein family, which are high ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MUC20
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	72
Uniprot ID	Q8N307
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...