


Product Datasheet

SENP2 antibody (orb66737)

Description	Rabbit polyclonal to SENP2.SENP2 is a protease that catalyzes two essential functions in the SUMO pa...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SENP2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	68
Uniprot ID	Q9HC62
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...