


Product Datasheet

SMARCAL1 antibody (orb66745)

Description	Rabbit polyclonal to SMARCAL1. SMARCAL1 is a member of the SWI/SNF family of proteins. It is an ATP-d...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SMARCAL1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	106
Uniprot ID	Q9NZC9
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...