


Product Datasheet

IL7Ra antibody (orb66759)

Description	Rabbit polyclonal to IL7Ra. Interleukin 7 receptor (IL7Ra) is the receptor for interleukine 7. The in...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse IL7Ra
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	52
Uniprot ID	P16872
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...