


Product Datasheet

CDH1 antibody (orb66765)

Description	Rabbit polyclonal to CDH1. Cadherin-1 (CDH1) is a single-pass membrane protein. Cadherins are calcium...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CDH1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	97
Uniprot ID	P12830
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...