

Product Datasheet

Integrin beta 1 antibody (orb66766)

Description	Rabbit polyclonal antibody to Integrin beta 1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse CD29
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	88
Uniprot ID	P09055
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

