


Product Datasheet

PRMT7 antibody (orb66768)

Description	Rabbit polyclonal to PRMT7. Protein arginine methyltransferase 7 (PRMT7) can catalyze the formation o...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PRMT7
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	78
Uniprot ID	Q9NVM4
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...