


Product Datasheet

SLC11A1 antibody (orb66774)

Description	Rabbit polyclonal to SLC11A1. SLC11A1 is a multi-pass membrane protein also known as Natural resistan...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SLC11A1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	60
Uniprot ID	P49279
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...