


Product Datasheet

IL25 antibody (orb66778)

Description	Rabbit polyclonal to IL25. Interleukin 25 (IL25) is a cytokine that induces activation of NF-kappa-B ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of human IL25
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	20
Uniprot ID	Q9H293
Dilution Range	WB:1:1000-1:5000, ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant hum...