


Product Datasheet

IL4 antibody (orb66783)

Description	Rabbit polyclonal antibody to IL4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of human IL4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	17
Uniprot ID	P05112
Dilution Range	WB:1:1000-1:5000,ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant hum...