


Product Datasheet

PSPC1 antibody (orb667833)

Description	Rabbit polyclonal antibody to PSPC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human PSPC1. The exact sequence is proprietary.
Target	PSPC1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8R326, Q8WXF1
Entrez	55269, 66645
Dilution Range	WB: 1-500-2000, IP: 1-20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of PSPC1 expression...